

Product datasheet for **HP214379**

CLSTN2 Human qPCR Primer Pair (NM_022131)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64084
Forward Sequence:	CACCATGTGGATGAAACACGGC
Reverse Sequence:	CTCCGCAAGAGAAAGACGAGG
ACCN:	NM_022131 , NM_022131.1 , NM_022131.2 , BC004871 , BC007943 , BC103496 , BC103506 , BC104485 , NM_022131.3
UniProt ID:	Q9H4D0
Synonyms:	ALC-GAMMA; alcagamma; CDHR13; CS2; CSTN2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)