

Product datasheet for **HP214387**

EPB41L4A Human qPCR Primer Pair (NM_022140)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64097
Forward Sequence:	CGAAGTTACGCCAGTATCGCA
Reverse Sequence:	CAGAACCTGAGCATCCGAAGA
ACCN:	<u>NM_022140</u> , <u>NM_022140.1</u> , <u>NM_022140.2</u> , <u>NM_022140.4</u> , <u>BC015441</u> , <u>BC031042</u> , <u>BC114632</u> , <u>BC114942</u> , <u>BC160044</u> , <u>BF431603</u> , <u>BX109898</u>
UniProt ID:	<u>Q9HCS5</u>
Synonyms:	EPB41L4; NBL4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)