

## Product datasheet for **HP214388**

### Parvin gamma (PARVG) Human qPCR Primer Pair (NM\_022141)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64098
Forward Sequence:	CTACCTCACTCCCAACTCTCCT
Reverse Sequence:	CTCAGTGTGCTCTTGGCATCCT
ACCN:	<a href="#">NM_022141</a> , <a href="#">NM_022141.1</a> , <a href="#">NM_022141.2</a> , <a href="#">NM_022141.3</a> , <a href="#">NM_022141.4</a> , <a href="#">NM_022141.5</a> , <a href="#">NM_022141.6</a> , <a href="#">BC034406</a> , <a href="#">BC034406.1</a> , <a href="#">BI818662</a> , <a href="#">BP299653</a> , <a href="#">BP299720</a> , <a href="#">BX648942</a>
UniProt ID:	<a href="#">Q9HBI0</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)