

## Product datasheet for **HP214401**

### DUS1L Human qPCR Primer Pair (NM\_022156)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64118
Forward Sequence:	GGAAGGAGAACCTGTACTGCGA
Reverse Sequence:	CCAAGTTCAGGTCAATGGCGTC
ACCN:	<u><a href="#">NM_022156</a></u> , <u><a href="#">NM_022156.1</a></u> , <u><a href="#">NM_022156.2</a></u> , <u><a href="#">NM_022156.3</a></u> , <u><a href="#">NM_022156.4</a></u> , <u><a href="#">BC064918</a></u> , <u><a href="#">BC064918.1</a></u> , <u><a href="#">BC003659</a></u> , <u><a href="#">BC017081</a></u> , <u><a href="#">BC062566</a></u> , <u><a href="#">BG709424</a></u> , <u><a href="#">BM907569</a></u> , <u><a href="#">NM_022156.5</a></u>
UniProt ID:	<u><a href="#">Q6P1R4</a></u>
Synonyms:	DUS1; PP3111
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)