

Product datasheet for **HP214412**

XYLT2 Human qPCR Primer Pair (NM_022167)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64132
Forward Sequence:	CCAACGACTTCAAGCCACAGGA
Reverse Sequence:	CTGCCATACAGGTGGAAGTCCA
ACCN:	NM_022167 , NM_022167.1 , NM_022167.2 , NM_022167.3 , BC052262 , BC052262.1 , BC017084 , BC036069 , BC156444 , BM726479 , BQ018370
UniProt ID:	Q9H1B5
Synonyms:	PXYLT2; SOS; XT-II; XT2; xylT-II
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)