

Product datasheet for **HP214424**

GLIPR2 Human qPCR Primer Pair (NM_022343)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	152007
Forward Sequence:	CCACTGAAGCTCTGCAAGAACC
Reverse Sequence:	AGGATGCCCATGCAAGGTTCTC
ACCN:	<u>NM_022343</u> , <u>NM_022343.1</u> , <u>NM_022343.2</u> , <u>NM_022343.3</u> , <u>BC017918</u> , <u>BF515272</u> , <u>BF516600</u> , <u>BP261391</u> , <u>NM_022343.4</u>
UniProt ID:	<u>Q9H4G4</u>
Synonyms:	C9orf19; GAPR-1; GAPR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)