

## Product datasheet for **HP214437**

### KCNK15 Human qPCR Primer Pair (NM\_022358)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	60598
Forward Sequence:	GGCTCCTTCTACTTCGCCATCA
Reverse Sequence:	CTGGAAAGTGACCAGCGTCAGC
ACCN:	<u><a href="#">NM_022358</a></u> , <u><a href="#">NM_022358.2</a></u> , <u><a href="#">NM_022358.3</a></u> , <u><a href="#">BC093874</a></u> , <u><a href="#">BC093874.1</a></u> , <u><a href="#">BC082981</a></u> , <u><a href="#">BC143282</a></u> , <u><a href="#">NM_022358.4</a></u>
UniProt ID:	<u><a href="#">Q9H427</a></u>
Synonyms:	dj781B1.1; K2p15.1; KCNK11; KCNK14; KT3.3; TASK-5; TASK5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)