

Product datasheet for **HP214442**

DNAJC1 Human qPCR Primer Pair (NM_022365)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64215
Forward Sequence:	GGACTCCAGGTCGATGGGAAAA
Reverse Sequence:	GGAGCAGGTCACCTGAATCCTTC
ACCN:	NM_022365 , NM_022365.1 , NM_022365.2 , NM_022365.3 , BC110894 , BC110894.1 , BC014436 , BC030955 , BC139752 , NM_022365.4
UniProt ID:	Q96KC8
Synonyms:	DNAJL1; ERdj1; HTJ1; MTJ1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)