

Product datasheet for **HP214445**

STRA6 Human qPCR Primer Pair (NM_022369)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64220
Forward Sequence:	GAATCGTGCTCTCCGAGGACAA
Reverse Sequence:	AGGACCAGGAAGGTGAGTAAGC
ACCN:	NM_022369 , NM_022369.1 , NM_022369.2 , NM_022369.3 , BC025256 , BC025256.1 , BC015881 , BX537413
UniProt ID:	Q9BX79
Synonyms:	MCOPCB8; MCOPS9; PP14296
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)