

## Product datasheet for **HP214447**

### TOR3A Human qPCR Primer Pair (NM\_022371)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64222
Forward Sequence:	GGTGGTCCTAAAGTTGCTCAAGG
Reverse Sequence:	CAGGTTTTCTTCACAAGACGGC
ACCN:	<u><a href="#">NM_022371</a></u> , <u><a href="#">NM_022371.1</a></u> , <u><a href="#">NM_022371.2</a></u> , <u><a href="#">NM_022371.3</a></u> , <u><a href="#">BC001085</a></u> , <u><a href="#">BC001085.1</a></u> , <u><a href="#">BC007571</a></u> , <u><a href="#">BC011746</a></u> , <u><a href="#">BC018292</a></u> , <u><a href="#">BC088368</a></u> , <u><a href="#">NM_022371.4</a></u>
UniProt ID:	<u><a href="#">Q9H497</a></u>
Synonyms:	ADIR; ADIR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)