

Product datasheet for **HP214448**

G protein beta subunit like (MLST8) Human qPCR Primer Pair (NM_022372)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64223
Forward Sequence:	CAGGTGAATGCCTTGAGGTCA
Reverse Sequence:	TTGTTGACGCCGTCGTAGCTGA
ACCN:	NM_022372 , NM_022372.1 , NM_022372.2 , NM_022372.3 , NM_022372.4 , BC088354 , BC088354.1 , BC001313 , BC017119 , BC020499 , BC052292 , BI761765 , BQ014640 , NM_022372.6
UniProt ID:	Q9BVC4
Synonyms:	GbetaL; GBL; LST8; POP3; WAT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)