

## Product datasheet for **HP214450**

### ARL6IP2 (ATL2) Human qPCR Primer Pair (NM\_022374)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64225
Forward Sequence:	GCAACGGTGTGGCTCTGAGCA
Reverse Sequence:	CATCGCAAGTCTCCATACTCTG
ACCN:	<a href="#">NM_022374</a> , <a href="#">NM_022374.1</a> , <a href="#">NM_022374.2</a> , <a href="#">NM_022374.3</a> , <a href="#">NM_022374.4</a> , <a href="#">BC039050</a> , <a href="#">BC053508</a> , <a href="#">BX507348</a>
UniProt ID:	<a href="#">Q8NHH9</a>
Synonyms:	aip-2; ARL3IP2; ARL6IP2; atlastin2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)