

Product datasheet for **HP214463**

Exportin 4 (XPO4) Human qPCR Primer Pair (NM_022459)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64328
Forward Sequence:	GCAGAGGACATTGATGAAGGCTC
Reverse Sequence:	AGAATCGCTGCTGAAGTGGCTG
ACCN:	<u>NM_022459</u> , <u>NM_022459.1</u> , <u>NM_022459.2</u> , <u>NM_022459.3</u> , <u>BC042504</u> , <u>BC050680</u> , <u>BC098402</u> , <u>BC113571</u> , <u>BC113573</u> , <u>BQ184183</u> , <u>BX647843</u>
UniProt ID:	<u>Q9C0E2</u>
Synonyms:	exp4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

