

## Product datasheet for **HP214464**

### HS1BP3 Human qPCR Primer Pair (NM\_022460)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64342
Forward Sequence:	GCCGTGTTCAATGAGATCCTGC
Reverse Sequence:	CAGAGGAATCTCTGCTGGTGAG
ACCN:	<a href="#">BC057389</a> , <a href="#">NM_022460</a> , <a href="#">NM_022460.1</a> , <a href="#">NM_022460.2</a> , <a href="#">NM_022460.3</a> , <a href="#">BC050636</a> , <a href="#">BC050636.1</a> , <a href="#">BC027947</a> , <a href="#">BC038847</a> , <a href="#">BG753647</a> , <a href="#">BX354596</a> , <a href="#">NM_022460.4</a>
UniProt ID:	<a href="#">Q53T59</a>
Synonyms:	ETM2; HS1-BP3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)