

## Product datasheet for **HP214468**

### ZNFN1A4 (IKZF4) Human qPCR Primer Pair (NM\_022465)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64375
Forward Sequence:	AGCCTCACCAAACGCAAGCGTT
Reverse Sequence:	AGGGAACTTCCAAAGCCAGGCT
ACCN:	<u>NM_022465, NM_022465.1, NM_022465.2, NM_022465.3, BC160169, BX647761, NM_022465.4</u>
UniProt ID:	<u>Q9H2S9</u>
Synonyms:	EOS; ZNFN1A4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)