

Product datasheet for **HP214491**

CIAO3 Human qPCR Primer Pair (NM_022493)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64428
Forward Sequence:	AGAGTTTGTGCGGCGATTCCGA
Reverse Sequence:	CGTGAGTCTTCTCGGCATAGCA
ACCN:	<u>NM_022493</u> , <u>NM_022493.1</u> , <u>NM_022493.2</u> , <u>BC030248</u> , <u>BC030248.1</u> , <u>BC129940</u> , <u>BU741568</u> , <u>NM_022493.3</u>
UniProt ID:	<u>Q9H6Q4</u>
Synonyms:	HPRN; IOP1; LET1L; NAR1; NARFL; PRN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)