

Product datasheet for **HP214499**

HLA-DRB3 Human qPCR Primer Pair (NM_022555)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3125
Forward Sequence:	TTCCAGACCCTGGTGATGCTAG
Reverse Sequence:	GACTCCACTCAGCATCTTGCTC
ACCN:	BC001023 , NM_022555 , NM_022555.1 , NM_022555.2 , NM_022555.3 , BC001023.2 , BC106057 , BC106057.1 , BC008987
UniProt ID:	P79483
Synonyms:	DRB3; HLA-DPB1; HLA-DR1B; HLA-DR3B; HLA-DRB3*
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)