

Product datasheet for **HP214503**

NDST4 Human qPCR Primer Pair (NM_022569)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64579
Forward Sequence:	TTGTGCAGAGCTGGACCAACCT
Reverse Sequence:	CTGGACCAGATGTCTTTGTGTCG
ACCN:	<u>NM_022569</u> , <u>NM_022569.1</u> , <u>NM_022569.2</u> , <u>BC012326</u> , <u>BC113075</u> , <u>BC156696</u> , <u>NM_022569.3</u>
UniProt ID:	<u>Q9H3R1</u>
Synonyms:	N-HSST; N-HSST 4; NDST-4; NHSST4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)