

Product datasheet for **HP214508**

VPS16 Human qPCR Primer Pair (NM_022575)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64601
Forward Sequence:	GTCTCTTACTCCGACATTGCTGC
Reverse Sequence:	GCTTAGTGCCAGTTTGCTCCTC
ACCN:	NM_022575 , NM_022575.1 , NM_022575.2 , NM_022575.3 , BC012422 , BC021291 , BC030571 , BC044662 , BC064406 , BC073959 , BC126431 , BM468030 , BU845142 , NM_022575.4
UniProt ID:	Q9H269
Synonyms:	DYT30; hVPS16
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)