

Product datasheet for **HP214538**

SFXN1 Human qPCR Primer Pair (NM_022754)

Product data:

| | |
|-------------------|---|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 94081 |
| Forward Sequence: | ACCAGTCCTTCAATGCCGTCGT |
| Reverse Sequence: | GAGTCCTAGAGCTGTTGCTACG |
| ACCN: | NM_022754 , NM_022754.2 , NM_022754.3 , NM_022754.4 , NM_022754.5 , NM_022754.6 , BC063241 , BC063241.1 , BC020517 , BP199230 , BU689165 , BX648188 , NM_022754.7 |
| UniProt ID: | Q9H9B4 |
| Synonyms: | SLC56A1; TCC |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)