

Product datasheet for **HP214554**

EPS8L2 Human qPCR Primer Pair (NM_022772)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64787
Forward Sequence:	CAAGGACCTGTTTGAGCAGAGG
Reverse Sequence:	GTGATGGCTTCGCTCTTGCCA
ACCN:	<u>NM_022772</u> , <u>NM_022772.1</u> , <u>NM_022772.2</u> , <u>NM_022772.3</u> , <u>BC080636</u> , <u>BC080636.1</u> , <u>BC002474</u> , <u>BC093878</u> , <u>BC101481</u> , <u>BC143242</u> , <u>NM_022772.4</u>
UniProt ID:	<u>Q9H6S3</u>
Synonyms:	DFNB106; EPS8R2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)