

## Product datasheet for **HP214570**

### LC3B (MAP1LC3B) Human qPCR Primer Pair (NM\_022818)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	81631
Forward Sequence:	GAGAAGCAGCTTCCTGTTCTGG
Reverse Sequence:	GTGTCCGTTACCAACAGGAAG
ACCN:	<u><a href="#">NM_022818</a></u> , <u><a href="#">NM_022818.2</a></u> , <u><a href="#">NM_022818.3</a></u> , <u><a href="#">NM_022818.4</a></u> , <u><a href="#">BC041874</a></u> , <u><a href="#">BC018634</a></u> , <u><a href="#">BC045759</a></u> , <u><a href="#">BC048009</a></u> , <u><a href="#">BC067797</a></u> , <u><a href="#">BG771621</a></u> , <u><a href="#">NM_022818.5</a></u>
UniProt ID:	<u><a href="#">Q9GZQ8</a></u>
Synonyms:	ATG8F; LC3B; MAP1A/1BLC3; MAP1LC3B-a
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)