

Product datasheet for **HP214576**

FBXL17 Human qPCR Primer Pair (NM_022824)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64839
Forward Sequence:	AGAGACTGTGGATGTCGGATGG
Reverse Sequence:	AGCTGTTCCACCGTCACTTCGT
ACCN:	<u>BC018548, NM_022824, NM_022824.1, BC126144, BC063316, BC126146, BM682767</u>
UniProt ID:	<u>Q9UF56</u>
Synonyms:	DKFZp434C1715; F-box and leucine-rich repeat protein 17; F-box only protein 13; Fbl17; Fbx13; FBXO13; MGC161422; MGC161424
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)