

Product datasheet for **HP214589**

METTL4 Human qPCR Primer Pair (NM_022840)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64863
Forward Sequence:	CTTGGTCTGTGGAGGTAGTTGC
Reverse Sequence:	CCAGTATAAGACCTTCGTAGGGC
ACCN:	<u>NM_022840</u> , <u>NM_022840.1</u> , <u>NM_022840.2</u> , <u>NM_022840.3</u> , <u>NM_022840.4</u> , <u>BC111020</u> , <u>BC056270</u> , <u>BC136766</u> , <u>BC136767</u> , <u>BI489580</u> , <u>BI771801</u> , <u>BM828675</u> , <u>BU741018</u> , <u>BX647373</u> , <u>NM_022840.5</u>
UniProt ID:	<u>Q8N3J2</u>
Synonyms:	HsT661
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)