

Product datasheet for **HP214591**

CDCP1 Human qPCR Primer Pair (NM_022842)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64866
Forward Sequence:	CAGGTGAAGCAGAACATCTCGG
Reverse Sequence:	GTCACCGTGAAAACGCCTTCCT
ACCN:	<u>NM_022842</u> , <u>NM_022842.1</u> , <u>NM_022842.2</u> , <u>NM_022842.3</u> , <u>NM_022842.4</u> , <u>BC021099</u> , <u>BC069254</u> , <u>BC148605</u> , <u>BC153169</u> , <u>NM_022842.5</u>
UniProt ID:	<u>Q9H5V8</u>
Synonyms:	CD318; SIMA135; TRASK
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)