

Product datasheet for **HP214617**

Nogo Receptor (RTN4R) Human qPCR Primer Pair (NM_023004)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	65078
Forward Sequence:	AGGCACTGCCTGATGACACCTT
Reverse Sequence:	TTCTGGTGCAGTAGGAGACGGT
ACCN:	NM_023004 , NM_023004.1 , NM_023004.2 , NM_023004.3 , NM_023004.4 , NM_023004.5 , BC011787 , BC011787.1 , BM762412 , NM_023004.6
UniProt ID:	Q9BZR6
Synonyms:	NGR; NOGOR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)