

Product datasheet for **HP214634**

FOXL2 Human qPCR Primer Pair (NM_023067)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	668
Forward Sequence:	CGGAGAAGAGGCTCACGCTGT
Reverse Sequence:	CTGAGGTTGTGGCGGATGCTAT
ACCN:	NM_023067 , NM_023067.1 , NM_023067.2 , NM_023067.3 , BC062549 , BC062549.1 , BC019896 , BF512713 , NM_023067.4
UniProt ID:	P58012
Synonyms:	BPES; BPES1; PFRK; PINTO; POF3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)