

## Product datasheet for **HP214666**

### CENPBD2P Human qPCR Primer Pair (NM\_023939)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	65996
Forward Sequence:	TGAGCCTCTTCCAGGAGTGGTT
Reverse Sequence:	AGCCAGACTCTGTCTCAACACC
ACCN:	<u><a href="#">NM_023939</a></u> , <u><a href="#">NM_023939.2</a></u> , <u><a href="#">BC000766</a></u> , <u><a href="#">BC001952</a></u> , <u><a href="#">BC025349</a></u> , <u><a href="#">BC043401</a></u> , <u><a href="#">BC065480</a></u> , <u><a href="#">BC067232</a></u> , <u><a href="#">BC098570</a></u> , <u><a href="#">BC119702</a></u> , <u><a href="#">BC119703</a></u> , <u><a href="#">BG214270</a></u>
Synonyms:	MGC111383, MGC149628
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)