

Product datasheet for **HP214672**

LYNX1 Human qPCR Primer Pair (NM_023946)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	66004
Forward Sequence:	ATCCAGATGCCTGAGGGACTCC
Reverse Sequence:	TGTGGCACATCTTGGTGACCAG
ACCN:	NM_023946.1 , NM_023946.2 , NM_023946.3 , BC014560 , BC047923 , BG701653 , BI603715 , BM172253 , BM671222 , BQ685632 , BU683767
UniProt ID:	P0DP58
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)