

Product datasheet for **HP214678**

IFNA13 (IFNA1) Human qPCR Primer Pair (NM_024013)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3439
Forward Sequence:	AGAAGGCTCCAGCCATCTCTGT
Reverse Sequence:	TGCTGGTAGAGTTCGGTGCAGA
ACCN:	NM_024013 , NM_024013.1 , BC074928 , BC074929 , BC112300 , BC112302 , NM_024013.3
UniProt ID:	P01562
Synonyms:	IFL; IFN; IFN-ALPHA; IFN-alphaD; IFNA13; IFNA@; IeIF D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)