

## Product datasheet for **HP214688**

### MRP63 (MRPL57) Human qPCR Primer Pair (NM\_024026)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	78988
Forward Sequence:	GGAGACGCAGAGTCTTGAGCAG
Reverse Sequence:	CGGTGCTCCCGATCCACTGA
ACCN:	<u><a href="#">NM_024026</a></u> , <u><a href="#">NM_024026.1</a></u> , <u><a href="#">NM_024026.2</a></u> , <u><a href="#">NM_024026.3</a></u> , <u><a href="#">NM_024026.4</a></u> , <u><a href="#">BC068492</a></u> , <u><a href="#">BC068492.1</a></u> , <u><a href="#">BC000002</a></u> , <u><a href="#">BC023616</a></u> , <u><a href="#">BG030730</a></u> , <u><a href="#">BG773624</a></u> , <u><a href="#">BG778041</a></u> , <u><a href="#">BI601227</a></u> , <u><a href="#">NM_024026.5</a></u>
UniProt ID:	<u><a href="#">Q9BQC6</a></u>
Synonyms:	bMRP63; MRP63
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)