

Product datasheet for **HP214700**

MIS12 Human qPCR Primer Pair (NM_024039)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79003
Forward Sequence:	AGTTACAGGAGAAGTACAAGACTG
Reverse Sequence:	CTAGTCCCATGATCTCTGCCAAC
ACCN:	NM_024039 , NM_024039.1 , BC000229 , BC000229.1
UniProt ID:	Q9H081
Synonyms:	2510025F08Rik; hMis12; KNTC2AP; MTW1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)