

Product datasheet for **HP214761**

RETREG2 Human qPCR Primer Pair (NM_024293)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79137
Forward Sequence:	GCTGTCAGATGAGGAGGCTTCT
Reverse Sequence:	TTCCTCTGGTTCACCTGGCAG
ACCN:	NM_024293 , NM_024293.1 , NM_024293.2 , NM_024293.4 , NM_024293.5 , BC002420 , BC064950 , BC115412 , BC141851 , BG700888 , BG723788 , BU615764 , NM_024293.6
UniProt ID:	Q8NC44
Synonyms:	C2orf17; FAM134A; MAG-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)