

Product datasheet for **HP214771**

APPD (PLEKHF1) Human qPCR Primer Pair (NM_024310)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79156
Forward Sequence:	TGGTGTATGGCAGCATCGTGCT
Reverse Sequence:	CGACACCACAAAGGACTTCTTGG
ACCN:	NM_024310 , NM_024310.1 , NM_024310.2 , NM_024310.3 , NM_024310.4 , BC002744 , BC002744.2 , BC018876 , NM_024310.5
UniProt ID:	Q96S99
Synonyms:	APPD; LAPF; PHAFIN1; ZFYVE15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)