

Product datasheet for **HP214773**

GNPTAB Human qPCR Primer Pair (NM_024312)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79158
Forward Sequence:	ATTCAGCCTGCGATTGGGATGG
Reverse Sequence:	TATTCCTCCACCAAAGTCCAGG
ACCN:	<u>NM_024312</u> , <u>NM_024312.1</u> , <u>NM_024312.2</u> , <u>NM_024312.3</u> , <u>NM_024312.4</u> , <u>BC131787</u> , <u>BC002779</u> , <u>BC042615</u> , <u>BC071687</u> , <u>BC094884</u> , <u>BP276948</u> , <u>BU754017</u> , <u>NM_024312.5</u>
UniProt ID:	<u>Q3T906</u>
Synonyms:	GNPTA; ICD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)