

## Product datasheet for **HP214810**

### ZKSCAN3 Human qPCR Primer Pair (NM\_024493)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	80317
Forward Sequence:	GCTAATGAAGGCTCTGCTCAAGC
Reverse Sequence:	GCCTAGAAGCTACCACTGCATC
ACCN:	<u><a href="#">NM_024493</a></u> , <u><a href="#">NM_024493.1</a></u> , <u><a href="#">NM_024493.2</a></u> , <u><a href="#">NM_024493.3</a></u> , <u><a href="#">BC006118</a></u> , <u><a href="#">BC006118.1</a></u> , <u><a href="#">BM762243</a></u> , <u><a href="#">BT007427</a></u>
UniProt ID:	<u><a href="#">Q9BRR0</a></u>
Synonyms:	dj874C20.1; dj874C20.1.; ZF47; zfp-47; Zfp47; ZFP306; ZNF306; ZNF309; ZSCAN13; ZSCAN35
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)