

Product datasheet for **HP214811**

IRF2BPL Human qPCR Primer Pair (NM_024496)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64207
Forward Sequence:	GAGGTTGCTTCAAGAAGGACC
Reverse Sequence:	GATGCACTGGAGTACACGTTGC
ACCN:	NM_024496 , NM_024496.1 , NM_024496.2 , NM_024496.3 , BC042460 , BC108292 , BC152431 , BI765830 , BM690017 , BM744266
UniProt ID:	Q9H1B7
Synonyms:	C14orf4; EAP1; NEDAMSS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)