

Product datasheet for **HP214824**

FYCO1 Human qPCR Primer Pair (NM_024513)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79443
Forward Sequence:	CCAGAAGATGCTTGCTGACCTG
Reverse Sequence:	AGCCATCTCTCCTCAGCACTGA
ACCN:	BC007218 , NM_024513 , NM_024513.1 , NM_024513.2 , NM_024513.3 , BC007218.1 , BC101468 , BC101470 , BC143368 , BP249484 , BX537863 , NM_024513.4
UniProt ID:	Q9BQS8
Synonyms:	CATC2; CTRCT18; RUFY3; ZFYVE7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)