

Product datasheet for **HP214840**

ERVMER34-1 Human qPCR Primer Pair (NM_024534)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	100288413
Forward Sequence:	GCTCACAGGTTCCAATCTGACC
Reverse Sequence:	GGTATCTGGTGAGGGAAGGTAC
ACCN:	<u>NM_024534</u> , <u>NM_024534.5</u> , <u>BC067526</u> , <u>BC067527</u> , <u>BC067528</u> , <u>BC067529</u> , <u>BC069998</u> , <u>BP871936</u>
UniProt ID:	<u>Q9H9K5</u>
Synonyms:	envMER34; HEMO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)