

Product datasheet for **HP214849**

C13orf7 (RNF219) Human qPCR Primer Pair (NM_024546)

Product data:

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 79596 |
| Forward Sequence: | TAGCCAGTGTCCAGCTTGCAGA |
| Reverse Sequence: | TCCGAAGATGCTTCTGACCGT |
| ACCN: | NM_024546 , NM_024546.1 , NM_024546.2 , BC028586 , BC136764 , BC136765 , BI550796 , BM674566 , NM_024546.4 |
| UniProt ID: | Q5W0B1 |
| Synonyms: | C13orf7; RNF219 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)