

Product datasheet for **HP214855**

RIC3 Human qPCR Primer Pair (NM_024557)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79608
Forward Sequence:	TACTGCCATGCCTGGAAACACC
Reverse Sequence:	CACCATTAGGTCCCACTCTGTTG
ACCN:	<u>NM_024557</u> , <u>NM_024557.1</u> , <u>NM_024557.2</u> , <u>NM_024557.3</u> , <u>NM_024557.4</u> , <u>NM_024557.5</u> , <u>BC022455</u> , <u>BC022455.1</u> , <u>BQ941375</u> , <u>BX647514</u>
UniProt ID:	<u>Q7Z5B4</u>
Synonyms:	AYST720; PRO1385
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)