

Product datasheet for **HP214867**

OGFRL1 Human qPCR Primer Pair (NM_024576)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79627
Forward Sequence:	TGGCAGGAAAGATTTTCAGCATCTG
Reverse Sequence:	AGTATTCTCCACAAGAGCTTCATG
ACCN:	NM_024576 , NM_024576.1 , NM_024576.2 , NM_024576.3 , NM_024576.4 , BC110999 , BC110999.1 , BC018688 , BX379107 , BX420584 , BX488500 , NM_024576.5
UniProt ID:	Q5TC84
Synonyms:	dj331H24.1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)