

Product datasheet for **HP214881**

SRD5A3 Human qPCR Primer Pair (NM_024592)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79644
Forward Sequence:	TCTTCCAGGACCTGATCCGCTA
Reverse Sequence:	CAGGAAGCCATTCCACAGCACT
ACCN:	NM_024592 , NM_024592.1 , NM_024592.2 , NM_024592.3 , NM_024592.4 , BC002480 , BC002480.2 , NM_024592.5
UniProt ID:	Q9H8P0
Synonyms:	CDG1P; CDG1Q; KRIZI; SRD5A2L; SRD5A2L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)