

Product datasheet for **HP214901**

TUT7 Human qPCR Primer Pair (NM_024617)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79670
Forward Sequence:	AGGCACCAAGAGAAACGCCGAT
Reverse Sequence:	CATAGAACCGCAGCAATTCCACC
ACCN:	<u>BC032456</u> , <u>NM_024617</u> , <u>NM_024617.1</u> , <u>NM_024617.2</u> , <u>NM_024617.3</u> , <u>BC032456.1</u> , <u>BC015897</u> , <u>BC016387</u> , <u>BC111003</u> , <u>BC144523</u> , <u>BU734789</u> , <u>BX641077</u> , <u>NM_024617.4</u>
UniProt ID:	<u>Q5VYS8</u>
Synonyms:	PAPD6; TENT3B; ZCCHC6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)