

Product datasheet for **HP214909**

B7H4 (VTCN1) Human qPCR Primer Pair (NM_024626)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79679
Forward Sequence:	CTCACAGATGCTGGCACCTACA
Reverse Sequence:	GCAAGGTCTCTGAGCTGGCATT
ACCN:	<u>NM_024626</u> , <u>NM_024626.1</u> , <u>NM_024626.2</u> , <u>NM_024626.3</u> , <u>BC074729</u> , <u>BC074729.2</u> , <u>BC065717</u> , <u>BX648021</u> , <u>NM_024626.4</u>
UniProt ID:	<u>Q7Z7D3</u>
Synonyms:	B7-H4; B7h.5; B7H4; B7S1; B7X; PRO1291; VTCN1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)