

Product datasheet for **HP214910**

RTL10 Human qPCR Primer Pair (NM_024627)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79680
Forward Sequence:	GTCCTGTGGAACCAGCCTCTTC
Reverse Sequence:	GGTTTTGGACCTCCTGGATGAG
ACCN:	NM_024627 , NM_024627.1 , NM_024627.2 , NM_024627.3 , NM_024627.4 , NM_024627.5 , BC011679 , BC011679.2 , BC008696 , BX640785 , BX640969 , BX640998 , NM_024627.6
UniProt ID:	Q7L3V2
Synonyms:	BOP; C22orf29
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)