

Product datasheet for **HP214920**

ZNF322A (ZNF322) Human qPCR Primer Pair (NM_024639)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79692
Forward Sequence:	CATCAGTGGTCTGCGTGTGAGA
Reverse Sequence:	GGCTGAACCTCTGGTGGAAGCTT
ACCN:	<u>NM_024639</u> , <u>NM_024639.2</u> , <u>NM_024639.3</u> , <u>NM_024639.4</u> , <u>BC050425</u> , <u>BC050425.1</u> , <u>BC074772</u> , <u>BC119669</u> , <u>BC139750</u> , <u>BX474628</u> , <u>NM_024639.5</u>
UniProt ID:	<u>Q6U7Q0</u>
Synonyms:	HCG12; ZNF322A; ZNF388; ZNF489
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)