

Product datasheet for **HP214946**

Lin28 (LIN28A) Human qPCR Primer Pair (NM_024674)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79727
Forward Sequence:	CCAGTGGATGTCTTTGTGCACC
Reverse Sequence:	GTGACACGGATGGATTCCAGAC
ACCN:	NM_024674 , NM_024674.1 , NM_024674.2 , NM_024674.3 , NM_024674.4 , NM_024674.5 , BC028566 , BC028566.2 , NM_024674.6
UniProt ID:	Q9H9Z2
Synonyms:	CSDD1; LIN-28; lin-28A; LIN28; ZCCHC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)