

Product datasheet for **HP214950**

NARS2 Human qPCR Primer Pair (NM_024678)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79731
Forward Sequence:	CTGTGGCAACATACCTGTCTTCG
Reverse Sequence:	CCAACTCCAGGAACCAGAAGATC
ACCN:	NM_024678 , NM_024678.1 , NM_024678.2 , NM_024678.3 , NM_024678.4 , NM_024678.5 , BC007800 , BC007800.2 , BM023108 , BP196275 , BP199309
UniProt ID:	Q96I59
Synonyms:	asnRS; DFNB94; SLM5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)