

Product datasheet for **HP214956**

TLL7 Human qPCR Primer Pair (NM_024686)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79739
Forward Sequence:	CTCTTACCAGGCTCCTCAGACT
Reverse Sequence:	GATCTGCTTTCGTA CTGAGCGA
ACCN:	NM_024686 , NM_024686.1 , NM_024686.2 , NM_024686.3 , NM_024686.4 , NM_024686.5 , BC028350 , BC032829 , BC042023 , BC048970 , BC060878 , BG718305 , NM_024686.6
UniProt ID:	Q6ZT98
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)